

BAI1 Rabbit pAb

Cat No: HR1APEA261

For research use only

Overview

Product Name	BAI1 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	130-170kDa
GeneID?Human?	575
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Brain specific angiogenesis inhibitor 1 antibody, GDAIF antibody, KIAA4089 antibody
Background	Brain-specific angiogenesis inhibitor 1 is a protein that in humans is encoded by the BAI1 gene.It is a member of the adhesion-GPCR family of receptors.