

Bestrophin-1 Rabbit pAb

Cat No: HR1APEA262

For research use only

Overview

Product Name	Bestrophin-1 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:50-100
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	67kDakDa
GeneID?Human?	7439
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BEST antibody, Bestrophin-1 antibody, Best disease antibody, VMD 2 antibody
Background	Bestrophins were first found by genetic linkage of human-Bestrophin-1 (hBest1) to a juvenile form of macular degeneration called Best vitelliform macular dystrophy (BVMD). Forms calcium-sensitive chloride channels.