

## Bombesin Receptor 2 Rabbit pAb Cat No: HR1APEA263

For research use only

## Overview

Product Name	Bombesin Receptor 2 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	
Clonality	Polycional
Concentration	1mg/ml
Observed band	50-70kDakDa
GenelD?Human?	2925
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BB2 antibody
Background	Bombesin receptor 2 (BB2) is a member of a family of receptors that binds the 14 amino acid peptide, bombesin.