

# Bombesin Receptor 2 Rabbit pAb

Cat No: HR1APEA263

For research use only

## Overview

|                       |                                                                                                                |
|-----------------------|----------------------------------------------------------------------------------------------------------------|
| Product Name          | Bombesin Receptor 2 Rabbit pAb                                                                                 |
| Source                | Rabbit                                                                                                         |
| Applications          | WB, IHC                                                                                                        |
| Species Reactivity    | Rat, Mouse                                                                                                     |
| Recommended Dilutions | WB 1:1,000-2,000 IHC 1:100-200                                                                                 |
| Immunogen             |                                                                                                                |
| Species               | Rabbit                                                                                                         |
| Storage               | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.        |
| Isotype               | IgG                                                                                                            |
| Clonality             | Polyclonal                                                                                                     |
| Concentration         | 1mg/ml                                                                                                         |
| Observed band         | 50-70kDa                                                                                                       |
| GeneID? Human?        | 2925                                                                                                           |
| Human Swiss-Prot No.  |                                                                                                                |
| Cellular localization |                                                                                                                |
| Alternative Names     | BB2 antibody                                                                                                   |
| Background            | Bombesin receptor 2 (BB2) is a member of a family of receptors that binds the 14 amino acid peptide, bombesin. |