

## Glutamate Receptor 1 Rabbit pAb Cat No: HR1APEA265

For research use only

## **Overview**

Product Name	Glutamate Receptor 1 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:50-100
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	100kDakDa
GeneID?Human?	2890
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GLUR1 antibody, GluA1 antibody, Glutamate receptor ionotropic antibody, GluRK1 antibody
Background	Glutamate receptor 1 is a protein that in humans is encoded by the GRIA1 gene. The receptor then desensitizes rapidly and enters a transient inactive state, characterized by the presence of bound agonist.