

CABP2 Rabbit pAb

Cat No: HR1APEA266

For research use only

Overview

Product Name	CABP2 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	37kDa
GeneID?Human?	51475
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DFNB93 antibody, CaBP2 antibody, Calcium-Binding Protein 2 antibody
Background	This gene belongs to a subfamily of calcium binding proteins that share similarity to calmodulin. Like calmodulin, these family members can likely stimulate calmodulin-dependent kinase II and the protein phosphatase calcineurin.