

Cav pan ?1 Rabbit pAb

Cat No: HR1APEA275

For research use only

Overview

Product Name	Cav pan ?1 Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	240kDa
GeneID?Human?	775
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CACNA1, CAC1C, CCHL1A1, DHPR alpha1,LQT8, TS, alpha-1 polypeptide
Background	Voltage-gated Ca ²⁺ channels (CaV), enable the passage of Ca ²⁺ ions in a voltage dependent manner. These heteromeric entities are formed in part by the pore-forming ?1 subunit which determines the biophysical and pharmacological properties of the channel.