

CCKBR Rabbit pAb

Cat No: HR1APEA276

For research use only

Overview

Product Name	CCKBR Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	50kDa
GeneID?Human?	887
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CCK B, CCKRB, Cholecystokinin 2 receptor,GASR, Gastrin receptor
Background	Receptor for gastrin and cholecystokinin. The CCK-B receptors occur throughout the central nervous system where they modulate anxiety, analgesia, arousal, and neuroleptic activity.