

For research use only

## Overview

| Product Name          | CLIC1 Rabbit pAb  |
|-----------------------|---|
| Source                | Rabbit  |
| Applications          | IHC   |
| Species Reactivity    | Human,Rat,Mouse   |
| Recommended Dilutions | IHC 1:100-200   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.                                 |
| Isotype               | IgG   |
| Clonality             | Polyclonal  |
| Concentration         | 1mg/ml  |
| Observed band         | 30kDakDa  |
| GeneID?Human?         | 1192  |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | Chloride channel ABP, NCC27,p64CLCP, RNCC,G6  |
| Background            | Membrane insertion seems to be redox-regulated and may occur only under oxydizing conditions. Involved in regulation of the cell cycle. |