

BTF3L4 Polyclonal Antibody

Cat No: HR1AP9078

For research use only

Overview

Product Name	BTF3L4 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	BTF3L4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BTF3L4; Transcription factor BTF3 homolog 4; Basic transcription factor 3-like 4
Background	similarity:Belongs to the NAC-beta family.,similarity:Contains 1 NAC-A/B (NAC-alpha/beta) domain.,