

CLIC4 Rabbit pAb
Cat No: HR1APEA278

For research use only

## **Overview**

Product Name	CLIC4 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	30kDakDa
GeneID?Human?	25932
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Chloride channel 4, FLJ38640,MC3S5,mtCLIC,p64H1,HUH1
Background	Membrane insertion seems to be redox-regulated and may occur only under oxydizing conditions. Involved in regulation of the cell cycle.