

CNTF Rabbit pAb
Cat No: HR1APEA279

For research use only

Overview

Product Name	CNTF Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG LIUEISU
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	30kDakDa
GeneID?Human?	1270
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Ciliary neurotrophic factor, HCNTF
Background	CNTF is a survival factor for various neuronal cell types. Seems to prevent the degeneration of motor axons after axotomy.