

## Connexin-26 Rabbit pAb Cat No: HR1APEA280

For research use only

## Overview

Product Name	Connexin-26 Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Mouse
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	26kDakDa
GeneID?Human?	2706
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GJB2,CX26,DFNA3,DFNB1, HID,NSRD1,PPK
Background	One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell.