

## Endothelin B Receptor Rabbit pAb

Cat No: HR1APEA282

For research use only

### Overview

Product Name	Endothelin B Receptor Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	50kDa
GeneID?Human?	1910
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ABCDS,Ednra,EDNRB,ET-B,ET-BR,HSCR2,WS4A
Background	N/A