

# EphA1 Rabbit pAb

Cat No: HR1APEA283

For research use only

## Overview

Product Name	EphA1 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	110kDa
GeneID?Human?	2041
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EPH Receptor A1, MGC163163,Tyrosin protein kinase receptor EPH
Background	The Eph receptors are the largest known family of receptor tyrosine kinases (RTKs). They can be divided into two groups based on sequence similarity and on their preference for a subset of ligands.