

GABA A Receptor α 2 Rabbit pAb

Cat No: HR1APEA285

For research use only

Overview

Product Name	GABA A Receptor α 2 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	51kDa
GeneID?Human?	282236
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GABA, GABRA2,GBRA2, GABA A receptor subunit alpha 2
Background	GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.