

GALR2 Rabbit pAb

Cat No: HR1APEA291

For research use only

Overview

Product Name	GALR2 Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	42kDa
GeneID?Human?	8811
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Galanin receptor 2,MGC125983
Background	GALR2 Receptor for the hormone galanin and for GALP. The activity of this receptor is mediated by G proteins that activate the phospholipase C/protein kinase C pathway (via G(q)) and that inhibit adenylyl cyclase (via G(i)).