

Ghrelin Receptor Rabbit pAb

Cat No: HR1APEA292

For research use only

Overview

Product Name	Ghrelin Receptor Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	41kDa
GeneID?Human?	2693
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GHRP, GHSR, GHRPR
Background	Receptor for ghrelin, coupled to G-alpha-11 proteins. Stimulates growth hormone secretion. Binds also other growth hormone releasing peptides (GHRP).