

## GLP1R Rabbit pAb

Cat No: HR1APEA293

For research use only

### Overview

Product Name	GLP1R Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	53kDa
GeneID?Human?	2740
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GLP, GLP1 receptor, Glucagon like peptide 1 receptor, MGC138331
Background	The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.