

KCNK10(TREK-2) Rabbit pAb

Cat No: HR1APEA294

For research use only

Overview

Product Name	KCNK10(TREK-2) Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	59kDakDa
GeneID?Human?	54207
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	K2p10.1,KCNKA,Potassium channel subfamily K member 10
Background	Outward rectifying potassium channel. Produces rapidly activating and non-inactivating outward rectifier K(+) currents. Activated by arachidonic acid and other naturally occurring unsaturated free fatty acids.