

## GRIK2(GluR6) Rabbit pAb

Cat No: HR1APEA297

For research use only

### Overview

Product Name	GRIK2(GluR6) Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	103kDa
GeneID?Human?	2898
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Excitatory amino acid receptor 4, GLR6, GLUK2, Glutamate receptor 6
Background	Ionotropic glutamate receptor. L-glutamate acts as an excitatory neurotransmitter at many synapses in the central nervous system.