

# KCNN2(SK2) Rabbit pAb

Cat No: HR1APEA300

For research use only

## Overview

Product Name	KCNN2(SK2) Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	70,26kDa
GeneID?Human?	3781
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Kca2.2,hsk2,skca2
Background	Forms a voltage-independent potassium channel activated by intracellular calcium. Activation is followed by membrane hyperpolarization.