

Cat No: HR1APEA303

For research use only

Overview

| Product Name | SLC12A4 Rabbit pAb |
|-----------------------|--|
| Source | Rabbit |
| Applications | WB, IHC |
| Species Reactivity | Human,Rat,Mouse |
| Recommended Dilutions | WB 1:1,000-2,000 IHC 1:100-200 |
| Immunogen | |
| Species | Rabbit |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Concentration | 1mg/ml |
| Observed band | 150kDakDa |
| GeneID?Human? | 6560 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | Electroneutral potassium chloride cotransporter 1,hKCC1,KCC1, s12A4,FLJ40489 |
| Background | Mediates electroneutral potassium-chloride cotransport when activated by cell swelling. May contribute to cell volume homeostasis in single cells. |