

SLC12A4 Rabbit pAb

Cat No: HR1APEA303

For research use only

Overview

Product Name	SLC12A4 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	150kDa
GeneID?Human?	6560
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Electroneutral potassium chloride cotransporter 1,hKCC1,KCC1, s12A4,FLJ40489
Background	Mediates electroneutral potassium-chloride cotransport when activated by cell swelling. May contribute to cell volume homeostasis in single cells.