

Kir4.1 Rabbit pAb

Cat No: HR1APEA306

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | Kir4.1 Rabbit pAb |
| Source | Rabbit |
| Applications | IHC |
| Species Reactivity | Human,Rat,Mouse |
| Recommended Dilutions | IHC 1:100-200 |
| Immunogen | |
| Species | Rabbit |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Concentration | 1mg/ml |
| Observed band | 42kDa |
| GeneID?Human? | 3766 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | BIRK10,KCNJ10,potassium channel |
| Background | N/A |