

KV?1 Rabbit pAb

Cat No: HR1APEA307

For research use only

Overview

Product Name	KV?1 Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	50kDa
GeneID?Human?	7881
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	hKvb3,KCNA1B,KV beta 1,Voltage-gated potassium channel beta 1 subunit
Background	N/A