

Kv1.3 Rabbit pAb

Cat No: HR1APEA308

For research use only

Overview

Product Name	Kv1.3 Rabbit pAb
Source	Rabbit
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	64kDa
GeneID?Human?	3738
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HGK5,KCNA3,MK3,Potassium channel 3, voltage gated potassium channel protein kv1.3
Background	N/A