

Kv11.1 Rabbit pAb

Cat No: HR1APEA313

For research use only

Overview

Product Name	Kv11.1 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Rat, Mouse
Recommended Dilutions	WB 1:1000-2000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	62kDa
GeneID? Human?	169522
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LQT2, HERG
Background	N/A