

Met Rabbit pAb Cat No: HR1APEA318

For research use only

Overview

Product Name	Met Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	130170kDakDa
GeneID?Human?	4233
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HGF, c-Met, tyrosine-protein kinase, oncogene met
Background	Met, a high affinity tyrosine kinase receptor for hepatocyte growth factor (HGF, also known as scatter factor) is a disulfide-linked heterodimer made of 45 kDa ?- and 145 kDa ?-subunits.