



NM23A Rabbit pAb

Cat No: HR1APEA322

For research use only

Overview

Product Name	NM23A Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	17kDa
GeneID?Human?	4830
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AWD,NME1,GAAD,NBS,NDKA,Nucleoside diphosphate kinase A antibody
Background	Members of the NDK/NME/NM23 kinase family inhibit metastasis in a variety of tumor cell types. All NDK/NME/NM23 proteins possess nucleoside diphosphatase kinase (NDK) activity and catalyze the phosphorylation of nucleoside diphosphate to the corresponding nucleoside triphosphate to regulate a diverse array of cellular events.