

AMPK γ 1 Rabbit pAb

Cat No: HR1APEA326

For research use only

Overview

Product Name	AMPK γ 1 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human, Mouse, Rat
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	38kDa
GeneID?Human?	5571
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AMPK gamma 1, AMPKg, PRKAG 1
Background	AMP-activated protein kinase (AMPK) is highly conserved from yeast to plants and animals and plays a key role in the regulation of energy homeostasis. AMPK is a heterotrimeric complex composed of a catalytic α subunit and regulatory β and γ subunits, each of which is encoded by two or three distinct genes (γ 1, 2; γ 1, 2; γ 1, 2, 3).