

## AMPK ?2 Rabbit pAb

Cat No: HR1APEA328

For research use only

## Overview

Product Name	AMPK ?2 Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	IHC 1:100
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	62kDakDa
GeneID?Human?	5563
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	5 AMP activated protein kinase catalytic subunit alpha-2, AAPK2, AMPK2, PRKAA2
Background	AMP-activated protein kinase (AMPK) is highly conserved from yeast to plants and animals and plays a key role in the regulation of energy homeostasis. AMPK is a heterotrimeric complex composed of a catalytic? subunit and regulatory? and? subunits, each of which is encoded by two or three distinct genes (?1, 2; ?1, 2; ?1, 2, 3).