

HP-1 ?/CBX1 Rabbit pAb

Cat No: HR1APEA333

For research use only

Overview

Product Name	HP-1 ?/CBX1 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	WB 1:1,000-5,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	25kDakDa
GeneID?Human?	10951
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CBX, Heterochromatin protein 1 beta, HP1 beta, p25 beta
Background	Heterochromatin protein 1 (HP1) is a family of heterochromatic adaptor molecules involved in both gene silencing and higher order chromatin structure. All three HP1 family members (?, ?, and ?) are primarily associated with centromeric heterochromatin.