

HSC70 Rabbit pAb Cat No: HR1APEA336

For research use only

Overview

Product Name	HSC70 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	70-72kDakDa
GeneID?Human?	3312
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Heat shock cognate 71 KDa protein antibody, HSC71 antibody, HSP73 antibody, HSPA8 antibody
Background	Heat shock 70 kDa protein 8 also known as heat shock cognate 71 kDa protein or Hsc70 or Hsp73 is a heat shock protein that in humans is encoded by the HSPA8 gene. As a member of the heat shock protein 70 family and a chaperone protein, it facilitates the proper folding of newly translated and misfolded proteins, as well as stabilize or degrade mutant proteins.