

HSP90? Rabbit pAb Cat No: HR1APEA338

For research use only

Overview

Product Name	HSP90? Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	95kDakDa
GeneID?Human?	3320
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EL52 antibody, Heat shock 86 KDa antibody, HSP86 antibody, HSP90N antibody, HSPC1 antibody
Background	Hsp90 (heat shock protein 90) is a chaperone protein that assists other proteins to fold properly, stabilizes proteins against heat stress, and aids in protein degradation. In mammalian cells, there are two or more genes encoding cytosolic Hsp90 homologues, with the human Hsp90? showing 85% sequence identity to Hsp90?.