

JAK1 Rabbit pAb

Cat No: HR1APEA340

For research use only

Overview

Product Name	JAK1 Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	130kDa
GeneID?Human?	3716
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Janus kinase 1,JTK3,Tyrosine protein kinase JAK1,JAK1A,JAK1B
Background	Members of the Janus family of tyrosine kinases (Jak1, Jak2, Jak3 and Tyk2) are activated by ligands binding to a number of associated cytokine receptors.