



JNK2 Rabbit pAb

Cat No: HR1APEA342

For research use only

Overview

Product Name	JNK2 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	WB 1:500-1,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	46,54kDa
GeneID?Human?	5601
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	c Jun kinase 2, JNK55, JNK2A, JNK2B, SAPK
Background	There are three SAPK/JNK genes each of which undergoes alternative splicing, resulting in numerous isoforms. SAPK/JNK, when active as a dimer, can translocate to the nucleus and regulate transcription through its effects on c-Jun, ATF-2, and other transcription factors.