

JNK3 Rabbit pAb

Cat No: HR1APEA343

For research use only

Overview

Product Name	JNK3 Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	54kDa
GeneID?Human?	5602
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	c Jun kinase 3, MAPK10, SAPK1b
Background	There are three SAPK/JNK genes each of which undergoes alternative splicing, resulting in numerous isoforms. SAPK/JNK, when active as a dimer, can translocate to the nucleus and regulate transcription through its effects on c-Jun, ATF-2, and other transcription factors.