

MEK2 Rabbit pAb

Cat No: HR1APEA344

For research use only

Overview

Product Name	MEK2 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	45kDa
GeneID?Human?	5605
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MAP2K2; MKK2; PRKMK2; Dual specificity mitogen-activated protein kinase kinase 2; MAP kinase kinase 2; MAPKK 2; ERK activator kinase 2; MAPK/ERK kinase 2
Background	MEK1 and MEK2, also called MAPK or Erk kinases, are dual-specificity protein kinases that function in a mitogen activated protein kinase cascade controlling cell growth and differentiation.