

ATG14L Rabbit pAb

Cat No: HR1APEA354

For research use only

Overview

Product Name	ATG14L Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	55kDa
GeneID?Human?	22863
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ATG14, KIAA0831,Autophagy-related protein 14-like protein antibody
Background	ATG14 is a autophagy-related protein. The serine/threonine kinase ULK1 phosphorylates Atg14 at Ser29 to promote autophagosome formation.