

FAK Rabbit pAb

Cat No: HR1APEA359

For research use only

Overview

Product Name	FAK Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	125kDa
GeneID?Human?	5747
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FADK, Focal adhesion kinase 1 antibody, Protein tyrosine kinase 2 antibody
Background	FAK is a widely expressed cytoplasmic protein tyrosine kinase involved in integrin-mediated signal transduction. It plays an important role in the control of several biological processes, including cell spreading, migration, and survival.