

Gab1 Rabbit pAb Cat No: HR1APEA360

For research use only

Overview

Product Name	Gab1 Rabbit pAb
Source	Rabbit
Applications	IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	76kDakDa
GeneID?Human?	2549
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GRB 2 associated binder 1 antibody
Background	The Gab1 protein is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis