

PGP9.5 (UCHL-1)Rabbit pAb

Cat No: HR1APEA365

For research use only

Overview

Product Name	PGP9.5 (UCHL-1)Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human, Mouse, Rat
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	25kDakDa
GeneID?Human?	7345
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Neuron cytoplasmic protein 9.5, PARK5, Ubiquitin thioesterase L1, Uch L1, UCHL1
Background	UCHL1, UCHL3, UCHL5/UCH37, and BRCA-1-associated protein-1 (BAP1) belong to the UCH family of DUBs, which all posses a conserved catalytic domain (UCH domain) of about 230 amino acids. UCHL1 is abundantly expressed in neuronal tissues and testes.