



Goat Anti-Rabbit IgG H&L Antibody (AP)

Cat No: HR1AP0010

For research use only

Overview

| | |
|-----------------------|-----------------------------------------------------------------------------|
| Product Name | Goat Anti-Rabbit IgG H&L Antibody (AP) |
| Source | Goat |
| Applications | WB |
| Species Reactivity | Rabbit |
| Recommended Dilutions | WB: 1:5000-1:50000 |
| Immunogen | |
| Species | Goat |
| Storage | Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Concentration | |
| Observed band | kDa |
| GenelD?Human? | |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | |