

Goat Anti-Rabbit IgG H&L Antibody (AP)

Cat No: HR1AP0010

For research use only

Overview

Product Name	Goat Anti-Rabbit IgG H&L Antibody (AP)
Source	Goat
Applications	WB
Species Reactivity	Rabbit
Recommended Dilutions	WB: 1:5000-1:50000
Immunogen	
Species	Goat
Storage	Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.
Isotype	IgG
Clonality	Polyclonal
Concentration	
Observed band	kDa
GeneID?Human?	
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	