

Rabbit IgG(Light Chain Specific) Mouse mAb Cat No: HR1AM2280

For research use only

Overview

Product Name	Rabbit IgG(Light Chain Specific) Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Rat
Recommended Dilutions	WB 1:5,000-20,000 IHC 1:200-500
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
lsotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	25kDakDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	The antibody reacts with the light chains on rabbit IgG and with those common to other rabbit immunoglobulins. The antibody does not react with the heavy chain of rabbit IgG. The antibody minimal cross-reaction with bovine, goat, Armenian hamster, horse, human, mouse, rat, and sheep immunoglobulins, but it may cross-react with IgG from other speciesThe antibody reacts with whole molecule rabbit IgG. The antibody may cross-react with immunoglobulins from other species.