

Anti-Baculovirus Envelope (gp64) antibody (monoclonal)

Cat No: HR6AM0085

For research use only

Overview

Product Name	Anti-Baculovirus Envelope (gp64) antibody (monoclonal)
Source	
Applications	IF, ELISA, FCM, ICC
Species Reactivity	Envelope Protein gp64 (baculovirus)
Recommended Dilutions	
Immunogen	
Species	
Storage	Store at -20°C or -80°C. Avoid repeated freeze.
Isotype	IgG2a
Clonality	Monoclonal
Concentration	
Observed band	kDa
GeneID?Human?	
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	