

Rabbit IgG(Fc specific, Heavy Chain) Mouse mAb

Cat No: HR1AM2279

For research use only

Overview

Product Name	Rabbit IgG(Fc specific,Heavy Chain) Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Rat
Recommended Dilutions	WB 1:5,000-20,000 IHC 1:200-500
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	55kDakDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	The antibody reacts with the Fc portion of rabbit IgG heavy chain but not with the Fab portion of rabbit immunoglobulins. The antibody may cross-react with immunoglobulins from other species.