



## TRMT11 Polyclonal Antibody

Cat No: HR1AP9475

For research use only

### Overview

Product Name	TRMT11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53kDa
GenID?Human?	TRMT11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TRMT11; C6orf75; MDS024; tRNA; guanine(10)-N2-methyltransferase homolog; tRNA guanosine-2'-O-methyltransferase TRM11 homolog
Background	function:Catalytic subunit of an S-adenosyl-L-methionine-dependent tRNA methyltransferase complex that mediates the methylation of the guanosine nucleotide at position 10 (m2G10) in tRNAs.,similarity:Belongs to the methyltransferase superfamily.,