




CGR19 Polyclonal Antibody

Cat No: HR1AP5207

For research use only

Overview

Product Name	CGR19 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GenelD?Human?	CGRRF1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CGRRF1; CGR19; RNF197; Cell growth regulator with RING finger domain protein 1; Cell growth regulatory gene 19 protein; RING finger protein 197
Background	function:Able to inhibit growth in several cell lines.,similarity:Contains 1 RING-type zinc finger.,



bioelsa